Stem Cell Work Station



Integrated Work Station for Stem Cell Isolation

Beauty Cell™ is an unique Multi-functional work station which combines a clean bench with a centrifuge and a shaking incubator. Its integrated design enables the operation of centrifuge and shaking incubator in clean bench which provides aseptic condition by air circulation through HEPA filter. It minimizes the risk of contamination during stem cell isolation or some biological application that using centrifuge and shaking incubator. Enjoy all-in-one station for the efficiency.









<UV Lamp>

<HEPA Filter>



<Inside Work Station>



<Centrifuge>

<Shaking Incubator>



<Interior Outlet>

<Cabinet for Storage>



<Pressure Gauze>



<Foot Switch>

Feature

#Clean Bench

It provides clean condition in work station when handling samples

- 1. **HEPA Filter** 99.99% efficiency for 0.3µm particles makes class 100 condition for safe cell handling. All airflow is through HEPA filer.
- 2. Equivlant Airflow as Biosafety Cabinet: 30% exhaust, 70% inflow
- 3. **Air Curtain**: Built-in air blower motor makes air curtain in front and prevent airflow from outside. It protects the sample from outside air.
- 4. **UV Lamp** for decontamination / **Fluorescent lamp** for lighting.
- 5. Storage Cabinet is located below work area for quick access on tools.
- 6. **Interior Power Outlet** for the use of an electric device in the station.
- 7. **Stainless steel work table** minimizes corrosion. Cleaning is also easy.
- 8. **Pressure Gauge** is to recognize the time for filter replacement.
- 9. Foot Switch is to operate power of centrifuge and shaking incubator

#Centrifuge


It is used for cell or blood isolation using centrifugal force.

- 1. Powerful but, stable centrifugation with quality inverter motor drive
- 2. **Teflon Coated Chamber** for easy cleaning and resistant on chemical
- 3. Microprocessor control, Unbalance Detection/ Buckets

#Shaking Incubator

It is used to mix and warm up the sample with enzyme or reagents.

- 1. Stable Shaking by BLDC Magnetic Induction motor
- 2. **Stable Temperature Control** by micro process
- 3. **Holder Platform** for 15, 50ml tubes and easy replaceable structure



Air Flow in the station

Adipose Stem Cell(SVF) Isolation Process




Liposuction



Preparation for SVF Isolation with Isolation kit



Enzymatic Digestion in shaking incubator



SVF Isolation by Centrifugation



Washing with Normal Saline and centrifugation



Final SVF Collection





Accessories

Swing Rotor Speed Table

Rotor	Bucket	VS-5500i		VS-5000i	
VR-50S	VB-04	5,500RPM	5,580xg	5,000RPM	4,612xg
VR-50S	VB-08	5,000RPM	4,612xg	4,500RPM	3,736xg
VR-50S	VB-32, VB-16	4,500RPM	3,736xg	4,000RPM	2,952xg
VR-50S	VC-40 carrier	3,500RPM	2,260xg	3,500PRM	2,260xg
VR-02M	Micro plate	2.500RPM	1.013xg	2,500RPM	1.013xp

ADSc(SVF) **Isolation Kit**

















Specification

Items	Unit	NB-803MS	NB-803MSF
Main Filter		HEPA filter 99,99% efficiency on particles of 0.3μm	HEPA filter 99,99% efficiency on particles of 0,3µm
Exhaust Filter		HEPA filter 99,99% efficiency on particles of 0,3μm	HEPA filter 99,99% efficiency on particles of 0,3µ
Sterilization lamp		U.V 20Wx2ea	U.V 20Wx2ea
CENTRIFUGE			
Max RPM		5,000rpm	5,500rpm
Max Force		4,612xg	5,580×g
Max Capacity	ml	480ml, 15mlx32Tubes(Swing rotor)	480ml, 15mlx32Tubes(Swing rotor)
Temp. range	°C		- 5℃~60℃
Main Controller		Digital PID controller	Digital PID controller
Drive system		Inverter Motor Drive	Inverter Motor Drive
Timer		99hrs 59min & hold run	99hrs 59min & hold run
SHAKING INCUBATOR	₹		
Temp range	°C	Ambient +5°C~60°C	Ambient +5°C∼60°C
Speed Range	rpm	30~300 rpm	30~300 rpm
Plate size	mm	250x310mm	250x310mm
Time Range		Continuous or up to 47h 59min	Continuous or up to 47h 59min
Display		LED Display	LED Display
Controller		Digital PID controller	Digital PID controller
DIMENSION			
Work space	mm	900x640x620(H)mm	1200x600x620(H)mm
Overall	mm	900x790x1920(H)mm	1200x825x1920(H)mm
Power	V/Hz	110/220V, 50/60Hz, 800W	110/220V. 50/60Hz. 1.2kW

© 2018 N-Biotek Inc. All rights reserved. Specifications, terms and appearance are subject to change.



Address: 402-803 Techno-Park, Pyeongcheon-ro 655, Wonmi-Gu Bucheon-Si, Republic of KOREA/ Zip Code 14502 E-mail: export@n-biotek.com Website: www.n-biotek.com